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
Description of material

Both **GumQuick** materials have been conceived on the basis of addition-vulcanising vinyl polysiloxanes for the production of flexible gingiva masks onto the master model.

Indication

GumQuick: Gingiva mask made in the dental laboratory for c & b.

GumQuick attains a final hardness of 60 Shore A and reproduces the natural proportions in combination with the opaque gingiva colour. Therefore unphysiological pressure onto the gingiva parts can already be avoided during the wax modellation.

GumQuick implant: Gingiva mask made in the dental laboratory for highly precise implant works.

GumQuick implant has been developed as special gingiva mask material for the implant technique. The high final hardness of 75 Shore A and the low flexibility of the material guarantee a safe and easy modellation of the secondary parts taking into consideration the natural course of the gingiva.

Technical Data according to DIN EN ISO 4823

	GumQuick	GumQuick implant
Processing time	1:30 Min.	1:30 Min.
Demoulding at 23 °C (incl. processing)	≥ 6:30 Min.	≥ 6:30 Min.
Tear-resistance	4,5 MPa	3,9 MPa
Elongation at break	175 %	55 %
Further tear-resistance	3,0 N/mm	2,0 N/mm
Recovery from deformation	> 99,5 %	> 99,5 %
Linear change of measures	< 0,5 %	< 0,5 %
Final hardness	60 Shore A	75 Shore A

The mixing and processing times refer to a room temperature of 23 °C and a relative air humidity of 50 %. Lower temperatures prolong, higher temperatures shorten these times.

This product was developed for dental use and has to be used according to the instructions. In the case of damages resulting from ulterior or improper use, the manufacturer assumes no liability. Furthermore the user is bound to check the product according

Processing instructions

GumQuick (indirect method)


A matrix of hard kneadable silicone material (e.g. **Laborsil pro**) is to be applied onto the dies of the unsawn die model which should be provided with the gingiva mask. After curing the kneadable silicone should be cut so that the top of the dies will appear (approx. 1 mm). The model can be sawn and gingiva parts ground as usual. Afterwards the correct position of the matrix has to be checked. In order to guarantee a faultless surface of the gingiva mask **Separator** should be applied onto the contact surface kneadable silicone / **GumQuick**.

GumQuick is to be applied directly from the cartridge onto the model and into the matrix. Now the matrix has to be positioned at the model so that the surplus of material can flow out through the holes at the top of the die. Ideally, the gingiva mask should be cured under pressure. After a period of approx. 6:30 minutes the matrix can be eliminated and the mask can be treated with scissors, scalpel or with rotating grinding instruments.

GumQuick implant (direct method)

At first **DuoSep** is to be applied onto the contact surface impression material / **GumQuick implant** in order to guarantee a faultless surface of the gingiva mask.

Attention: In order to avoid a damage of the gingiva mask surface when applying **GumQuick implant** in conjunction with polyether impressions **DuoSep** is to be used! Before the production of the gingiva mask the polyether impressions should ventilate for at least 3 hours. Afterwards **GumQuick implant** is to be applied directly from the cartridge into the impression in the area of the impression posts. Thanks to the ideal flowing characteristics of **GumQuick implant** the abutments will be surrounded

safely. After a period of approx. 6:30 minutes **GumQuick implant** is cured and the impression can be filled with plaster. The gingiva mask can easily be treated with scissors, scalpels or with rotating grinding instruments.

Important advice

Due to technical filling reasons an unequal transportation of the components can occur when using a new cartridge. For this reason the first transported quantity – approx. size of a pea – should not be used. Never try to squeeze the material through cured mixing canulas. This may lead to leakages in the cartridge system and therefore to the emission of unmixed and uncured single components. Before each application the cartridge has to be checked for consistent transportation. The vulcanisation can be disturbed by the use of latex gloves, hand creams, cleaning agents etc. The product is exclusively designed for the a. m. field of application and should only be used by qualified and instructed persons.

Keep away from children!

Only for the a. m. field of application!

Storage conditions

At room temperature, i. e. 18 °C – 28 °C.

Durability: 2 years in the closed original container.

GumQuick / GumQuick implant is not to be used when expiry date has passed.

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