GB - Henry Schein® Vinylpolysiloxane (VPS Hydro) Impression Material Light Body, Fast Set and Light Body, Regular Set, Regular Body Instructions for Use**Description:**

ISO 4823 Type 3 Light-bodied Consistency - ADA 19 Type 1 Low Viscosity
ISO 4823 Type 2 Medium-bodied Consistency - ADA 19 Type 1 Medium Viscosity

Henry Schein® Vinylpolysiloxane (VPS Hydro) is a hydrophilic addition reaction silicone impression material dispensed in an automix cartridge delivery system. Henry Schein® VPS Hydro Impression Material is thixotropic - it resists slumping, yet flows under pressure. Henry Schein® VPS Hydro Light Body is available in two setting types: Fast Set and Regular Set. Light Body, Fast Set is designed primarily for use with 1 or 2 unit impressions where a faster setting time is required.

Indications: Correction impression material for the two-step (putty-wash) impression technique. Use with dispenser to extrude. Use Henry Schein® Intra-Oral Tips for extruding material directly into mouth.

Impression Procedure: Inject automix material evenly around prepared teeth and/or directly onto Putty Impression. Seat tray and maintain in position for approximately 2,30 minutes (Fast Set) or 3,30 minutes (Regular Set). **Note:** Set time is based on 23°C. Higher temperature decreases set time, lower temperature increases set time.

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Electroplating: Impressions may be copper or silverplated.

Storage: Do not store above 27°C (80°F). Shelf life 3 years when stored at 23°C (73°F).

General Notes: Do not use with condensation reaction silicone impression materials.

TECHNICAL DATA

	REGULAR SET	FAST SET
Mixing Ratio Base : Catalyzt	1:1	1:1
Mixing Time	Auto Mixing	Auto Mixing
Total Working Time* (includes Mixing Time)	2 min	1:30 min
Time in mouth* (minimum)	3:30 min	2:30 min
Setting Time* (includes Mixing Time)	5:30 min	4 min
Linear Dimensional Change after 24 Hours	<0,2%	<0,2%
Maximum Strain in Compression	3-5 %	3-5 %
Recovery from Deformation	>99,5%	>99,5%
Hardness "Shore A" *23°C/73°F	45 ± 2	45 ± 2

Disinfection: Impression may be disinfected or sterilized with a disinfectant for impressions.

Compatibility with Gypsum: This product is compatible with all well known gypsum materials (type 3 Gibraltar Stone - Henry Schein; type 4 Gibraltar Die Stone - Henry Schein).

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Hydrophilic/Hidrófilo/Hydrophile/Hydrophil

HENRY SCHEIN®

Light Body, Fast Set
Light Body, Regular Set
Regular Body

VPS Hydro

Vinylpolysiloxane Impression Material
Light Body, Fast Set, Light Body, Regular Set & Regular Body
Instruction for use

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GB - Henry Schein® Vinylpolisiloxane (VPS) Putty Instructions for Use

Description:
ISO 4823 Type 0 Putty Consistency - ADA 19 Type 1 Very High Viscosity Henry Schein® Vinylpolisiloxane (VPS) Hydro putty is a hydrophilic addition reaction silicone impression material.

Indications: Preliminary impression for the two-step (putty-wash) technique in combination with Henry Schein® VPS Hydro Light Body, Regular Set/Fast Set or Regular Body syringe materials.

Trays and Adhesives: Perforated or non-perforated metal or plastic trays may be used with Henry Schein® VPS Hydro putty. To improve adhesion of Putty to the tray, Henry Schein® VPS Adhesive should be applied in a thin, even layer. Allow adhesive to dry for 3 minutes before loading the impression tray.

Dispensing and Mixing: Dispense equal volumes or weights of catalyst and base, using the colored spoon for the base and the white spoon for the catalyst. Do not interchange spoons. Knead with fingertips until uniform and streak-free (approximately 30 seconds). **Note:** Containers should be closed immediately after use, taking care not to interchange lids.

Impression Procedure: Load properly mixed material into the prepared tray. Place a polyethylene spacer sheet over putty and seat tray in the mouth not longer than 1 minutes after start of mix. Maintain position for approximately 2 minutes (material should have a rubbery consistency). **Note:** Higher temperature decreases set time; conversely, lower temperature increases set time.

Casting: Store impression at room temperature. Impression may be poured 30 minutes after removal from mouth.

Electroplating: Impressions may be copper or silver plated.

Storage: Do not store above 27°C (80°F). Shelf life is 3 years when stored at 23°C (73°F).

- General**
- Do not interchange catalyst and base covers. Do not interchange measuring spoons.
 - Certain gloves will inhibit the setting of putty. Test gloves to confirm proper setting times. Vinyl gloves or vinyl overliners are recommended. DO NOT USE LATEX GLOVES.
 - Do not use with condensation reaction silicone impression materials.

TECHNICAL DATA	REGULAR SET	FAST SET
Mixing Ratio Base : Catalyst	1:1	1:1
Mixing Time	30 s	30 s
Total Working Time* (includes Mixing Time)	2 min	1:30 min
Time in mouth* (minimum)	3:30 min	2:30 min
Setting Time* (includes Mixing Time)	5:30 min	4:00 min
Linear Dimensional Change after 24 Hours	<0.2%	<0.2%
Maximum Strain in Compression	1-3%	1-3%
Recovery from Deformation	>99%	>99%
Hardness "Shore A"	60 ± 2	60 ± 2

Disinfection: Impressions may be disinfected or sterilized with a disinfectant for impressions.

Compatibility with Gypsums: This product is compatible with all well-known gypsum materials (type 3 Gibraltar Stone - Henry Schein; type 4 Gibraltar Die Stone - Henry Schein).

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Henry Schein U.K. Holdings Ltd.