## ACRYL COLD (GB)

**Indications:** Adding the saddles to chrome cobalt denture bases. Repairing, extending and relining dentures using liquid resin.

**Properties:** Color stability, high mechanical values

**Advantages:** Easily processed.

**Mixing ratio:** 10 parts by weight powder, 7 parts by weight liquid.

**Processing:** Apply separating agent to the plaster. Add the powder to the liquid, allow 15 seconds for it to soak in completely and spatulate the Acryl Cold thoroughly in a mixing bowl. In its fluid state, Acryl Cold can be poured into, e. g. an overcast. Once the plastic phase sets in, the Acryl Cold will no longer flow out of the overcast and can be shaped. To prevent acrylic teeth breaking off the denture, the teeth must be roughened, as it is standard procedure for cold-curing resins, provided with undercuts and wetted with monomer. After a maximum of 6 minutes (at a room temperature of 23°C/73.4°F), Acryl Cold must be placed in a pressurized polymerization system for polymerization.

**Polymerization:** ACRYL COLD should be cured in a pressurized polymerization system for 15 minutes, at a pressure of 2-5 bar and a water temperature of 45°C/113°F.

**Repairs - Please note:** Repairs can be carried out with Acryl Cold. Prior to applying the resin dough, the fractured surfaces must be roughened, cleaned and wetted with monomer. The resin should be cured in a pressurized polymerization system for 15 minutes, at a pressure of 2-5 bar and a water temperature of 45°C/113°F.

**Storage:** Keep at room temperature - approximately 21°C/69.8°F.

**Important:** Should the room temperature deviate from 23°C/73.4°F: Higher temperatures reduce and lower temperatures prolong the working times.

**Colors:** Pink-transparent.

**Presentation:** 500g and 1000g Powder or 500ml and 1000ml Liquid.

**Danger:** Avoid prolonged skin contact with unpolymerized material and inhalation of monomer vapours. Once the container has been opened: Do not use after the expiry date.

**Disposal:** The powder can be disposed of with household waste. Dispose of the liquid in a special refuse tip.

**Side-effects:** With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects - even if it is doubtful that the side-effect has been caused by our product - please kindly contact us.

**Contra-indications/interactions:** If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth before using the product.

**Note:** Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

PROBLEM	TROUBLESHOOTING CAUSE	REMEDY
Acrylic does not fully cure	Incorrect mixing ratio	Adhere to the instructions for use Adhere to the mixing ratio
Acrylic is whitish on the fitting surface	Incorrect or poor application of separating agent	Apply ample separating agent Check separating agent if necessary Soak the model for approx. 10 min. in hand-hot water
Acrylic form bubbles on the fitting surface	Model not soaked	Soak the model for approx. 10 min. in hand-hot water
Acrylic surface white/porous	Model not soaked, recommended time for injecting acrylic exceeded or delay in inserting into the pressure vessel	Adhere to the processing times
Colour differences with repairs	Incorrect mixing ratio	Adhere to the mixing ratio
Acrylic breaks when lifted	Poorly separated	Check the separating agent
Acrylic/denture has whitish streaks	Incorrectly mixed/too much powder used	Thoroughly spatulate the acrylic dough/adhere to the mixing ratio

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Henry Schein U.K. Holdings Ltd.