## INSTRUCTIONS FOR USE

### Henry Schein Resin Cement SE

HS-Resin Cement SE is permanent resin cement does not require etching, priming or bonding. Self / light cured. The "Automix" system provides a rapid, consistent and homogeneous mixture which can be dispensed directly into the resation or tooth cavity.

**INGREDIENTS:**

***Base & Catalyst:***

- 4-Methacryloxyethyltrimellitic anhydride .
- Multifunctional monomers.
- Urethan dimethacrylate.
- Photoinitators.
- Polymerization activator.
- Glass filler.
- Fumed Silica.

**INSTRUCTIONS:**

- Remove the temporary restoration. Remark: Please only use temporary cements without eugenol. Eugenol-containing materials products can inhibit the polymerization process
- Rinse carefully
- Dry the preparation with air (do not dry up).
- Try in the restoration to ensure proper fit.
- Bleed the syringe before initial use. Place mixing tip on the Automix syringe.
- Dispense the material as follows: Inlay/Onlay: Mount the supplied intra oral tip onto the mixing tip and dispense the cement directly covering all surfaces. Put the restoration gently onto the preparation allowing the material to flow from all sides.
- Crown: Dispense the cement directly into the restoration covering all surfaces. Place the restoration smoothly onto the preparation allowing the cement to flow from all margins.
- Post: Mount the supplied intra oral tip onto the mixing tip. Put the cement into the canal. Place the post and move slightly avoiding trapped air. Help the cement to flow from the canal space. When the restoration is properly placed, please remove all extra cement. Excess cement is best to be removed when it is still in gel stage with a scaler / explorer. Gel consistency can be achieved by curing excess with a light for 2 seconds only or allowing to self-cure for 1.5 - 2 minutes after application. To achieve the best results, leave it undisturbed for 1.5 minutes before light curing.

**TIMES:**

- Working time: 2 minutes.
- Oral setting time: 3 minutes.
- Light-curing time: 20 seconds (see manufacturer's recommendation).
- After removal of excess cement, light cure all surfaces including margins for 20 seconds.

Use finishing strips for interproximal areas. Polish margins with polishing discs, cups, or points. Check occlusion.

**Remark:** Automix tips are designed for single patient use only.

**STORAGE:** Store below 25 °C (77 °F). Keep away from direct sunlight. For optimum freshness, keep refrigerated. Do not freeze. If refrigerated, allow the syringe to reach room temperature. HS-Resin Cement SE is not recommended for use beyond it's stated expiration date marked on the packaging. Shelf life is 2 years after the date of manufacturing.

**CAUTION!** Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact.

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HENRY SCHEIN SERVICES GmbH