Consult, consult accompanying documents

Do not reuse

Medical Device

Item Code	Product Description
9011255	ZIRLUX ACETAL MILLING DISC A1 98.5 x 15mm
9011260	ZIRLUX ACETAL MILLING DISC A1 98.5 x 20mm
9011266	ZIRLUX ACETAL MILLING DISC A1 98.5 x 25mm
9011256	ZIRLUX ACETAL MILLING DISC A2 98.5 x 15mm
9011261	ZIRLUX ACETAL MILLING DISC A2 98.5 x 20mm
9011267	ZIRLUX ACETAL MILLING DISC A2 98.5 x 25mm
9012701	ZIRLUX ACETAL MILLING DISC A3 98.5 x 15mm
9011257	ZIRLUX ACETAL MILLING DISC A3.5 98.5 x 15mm
9012702	ZIRLUX ACETAL MILLING DISC A3 98.5 x 20mm
9011262	ZIRLUX ACETAL MILLING DISC A3.5 98.5 x 20mm
9012703	ZIRLUX ACETAL MILLING DISC A3 98.5 x 25mm
9011268	ZIRLUX ACETAL MILLING DISC A3.5 98.5 x 25mm
9011258	ZIRLUX ACETAL MILLING DISC B1 98.5 x 15mm
9011263	ZIRLUX ACETAL MILLING DISC B1 98.5 x 20mm
9011269	ZIRLUX ACETAL MILLING DISC B1 98.5 x 25mm
9012704	ZIRLUX ACETAL MILLING DISC C3 98.5 x 15mm
9012705	ZIRLUX ACETAL MILLING DISC C3 98.5 x 20mm
9012706	ZIRLUX ACETAL MILLING DISC C3 98.5 x 25mm
9011259	ZIRLUX ACETAL MILLING DISC G2 (BLEACHED) 98.5 x 15mm
9011264	ZIRLUX ACETAL MILLING DISC G2 (BLEACHED) 98.5 x 20mm
9011270	ZIRLUX ACETAL MILLING DISC G2 (BLEACHED) 98.5 x 25mm
9012436	ZIRLUX ACETAL MILLING DISC P1 98.5 x 15mm
9012437	ZIRLUX ACETAL MILLING DISC P1 98.5 x 20mm
9012438	ZIRLUX ACETAL MILLING DISC P1 98.5 x 25mm
9012490	ZIRLUX ACETAL MILLING DISC A1 71 x 20mm
9012491	ZIRLUX ACETAL MILLING DISC B1 71 x 20mm B1
9012492	ZIRLUX ACETAL MILLING DISC G2 (BLEACHED) 71 x 20 mm

Zirlux Acetal Instructions for Use

SECTION 1: SYSTEM OVERVIEW

Zirlux Acetal is a pure, highly crystalline acetal copolymer resin possessing high tensile and flexural strength, fatigue resistance, and hardness. These properties, when combined with its very low moisture absorption, make the acetal material ideally suited for milling a variety of removable dental appliances, including partial denture frameworks, provisional crowns and bridges, and bite splints.

SECTION 2: PRODUCT COMPOSITION

Ingredient	Approx % weight
Polyoxymethylene Copolymer	99 to 100
Pigments	< 1.0

SECTION 3: TECHNICAL DATA

Property	Method	Unit	Value
Density	ASTM D792	g/cm ³	1.41
Tensile Strength (New)	ASTM D638	MPa	61.00
Tensile Modulus	ASTM D638	MPa	2820.00
Flexural Strength (New)	ASTM D790	MPa	90.00
Flexural Modulus	ASTM D790	MPa	2620.00
Moisture Absorption	ASTM D570	wt%	0.22
Izod Impact, Notched	ASTM D256	ft-lb/in (J/m)	1.00

SECTION 4: INDICATIONS

Removable prosthetic dental appliances
Partial denture clasp
Partial denture framework
Provisional crown & bridge
Bite splint

SECTION 5: CONTRAINDICATIONS

Contraindications are all applications not included in SECTION 4 above.

SECTION 6: PREPARATION GUIDELINES

General preparation guidelines for all removable appliances should be followed.

SECTION 7: DESIGN OF RESTORATION

Partial denture framework (minimum thickness)

- Major connector: 0.7 mm
- Saddle area: 0.5 mm
- Clasp on anterior tooth
- 2.0 mm where clasp joins framework, thinning to 0.7 mm at tip

Clasp on posterior tooth

- 2.5 mm where clasp joins framework, thinning to 0.7 mm at tip

Temporary Crown & Bridge

- 0.5 mm minimum thickness
- 0.03 mm cement gap

SECTION 8: MILLING

1. Mount the disc according to the CAM-System instructions. Set the enlargement-factor into the unit at a ratio of 1:1.
2. Start the milling unit.
3. Select appropriate milling strategy depending on the appliance to be fabricated.
4. Once the milling operation has completed, remove the disc containing the milled restoration(s).

SECTION 9: POST MILLING

1. Remove the restoration(s) from the disc using a separating disc or suitable acrylic bur.
2. To polish the restoration, load a new polishing wheel or buff with polish suitable for thermoplastic materials and buff all areas of the appliance until all scratches are removed and high shine is achieved. Notice: All Zirlux Acetal discs are solely for use by or on the order of a dental professional. They are not for use by the general public or over-the-counter.

SECTION 10: HANDLING & STORAGE

Please refer to MSDS.

SECTION 11: REGULATORY

ISO Specifications 9001 and 13485.

According to the EU Medical Devices Regulation, users / patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred.

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