BACS125
BACS126

English / Deutsch / Français / Español / Italiano / Nederlands

Operating Instructions
Gebrauchsanweisung
Notice d'utilisation

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Operating Instructions

Table of Contents

1	Before you begin	5
1.1	Structure of the document	6
1.2	Service life	9
2	Safety instructions	10
3	Technical description	14
3.1	Task	14
3.2	Ultrasonic handpiece design.....	14
3.3	Technical Data.....	15

4	Pre-reprocessing.....	18
4.1	Initial start-up and longer breaks in use.....	18
4.2	Prior to starting the work day.....	18
4.3	Before each patient.....	18
5	Operation	20
5.1	Connecting the ultrasonic handpiece to the supply hose	21
5.2	Detaching the ultrasonic handpiece from the supply hose	22
5.3	Inserting and removing the ultrasonic tip.....	23
5.4	Adjusting the cooling spray.....	25
5.5	Treatment	26
6	Post-reprocessing	29
6.1	After each treatment session.....	29

7	Reprocessing	31
7.1	Conduct pre-disinfection	31
7.2	Automated cleaning and disinfection	32
7.3	Manual cleaning and disinfection.....	34
7.4	Sterilizing	37
8	Maintenance	40
8.1	Maintaining the water channel	40
8.2	Cleaning the light guide surfaces.....	41
8.3	Replacing the spray water cartridge	41
8.4	Replacing the halogen lamp/LED	43

1 Before you begin ...

BACS125 / BACS126 complies with the regulations taking into account the state of the art.

1. Read the operating instructions prior to using BACS125 / BACS126.
2. Use BACS125 / BACS126 only for the applications described in the operating instructions.
3. Observe the hygiene standards, occupational safety regulations, and accident prevention measures applicable to the BACS125 / BACS126 in your country.

BACS125 / BACS126 is used to treat dental disease. The treatment is for the following uses:

- Supragingival removal of plaque and tartar (scaling)
- Subgingival removal of plaque and calculus (periodontology)

Intended use

- Cleaning implants as well as metallic or ceramic restorations (implant cleaning tip)
- Irrigation of prepared root canals (endo),
- liquefying of thixotropic composites (CEM),
- Micropreparation
- Retrograde root treatment (endo)

none

This product is intended only for use by trained dental personnel in dental practices and laboratories.

Contraindications

Target group

1.1 Structure of the document

1.1.1 Labeling of information

➤ To prevent injuries, please observe the warnings.

Warnings are labeled as follows:

Warnings

Instructions for use

⚠ DANGER! indicates a danger **leading** to death or serious injury if not avoided.

⚠ WARNING! indicates a danger that **may lead** to death or serious injury if not avoided.

⚠ CAUTION! indicates a danger that **may lead** to injury if not avoided.

➤ To prevent material damage and additional expenses, please observe all instructions for use.

Instructions for use are labeled as follows:

NOTE! indicates measures for the prevention of material damage.

IMPORTANT: indicates information on the avoidance of additional expenses and other important information.

Tip: indicates information for facilitating work.

English

1.1.2 Formats and symbols

The formats and symbols used in this document have the following meaning:

<ul style="list-style-type: none"> ✓ Prerequisite 1. First action step 2. Second action step or ➤ Alternative action ↪ Result ➤ Individual action step 	Requests you to do something.
Use of formats and symbols [→ 8].	Identifies a reference to another text passage and specifies its page number.
• List	Identifies a list.

1.2 Service life

When used as intended:

- Non-moving parts have a typical service life of 5 years
- Moving parts have a typical service life of 3 years

No warranty claim can be inferred here, as wear may occur earlier or later than indicated above depending on use, frequency of sterilization, and frequency of maintenance.

Obligations of the user

2 Safety instructions

- Use only fault-free materials that do **not** deviate from the specified data [→ 15].
- Protect yourself, patients, and others from danger. To do this, observe the safety information.
- Comply with the Intended use of the equipment.
- You should always keep these operating instructions within reach for further reference.
- Use BACS125 / BACS126 only with devices approved by the manufacturer.
- Use only the torque wrench included with each tip.
- Dispose of the torque wrench with the tip after using them for 12 months.

Preventing the spread of infections and cross contamination

Prevent the spread of infections and cross contamination between patients, users, and third parties. Sterilize equipment after each patient.

Take the appropriate hygiene measures, e.g. wear protective gloves.

Prevention of eye damage

The LED is in risk class 2 according to the IEC 62471:2006 standard. The LED emits optical radiation that is potentially hazardous and may be harmful to the eyes! Potential damage to the retina from the blue light emission. Do not stare at the LED for longer periods of time while in operation.

Improper use

Improper use may cause damage to the surface of the patient's teeth.

Malfunction or damage

Discontinue use immediately in case of malfunction, unusual or different sounds or damage. Damaged ultrasonic handpieces and tips may cause injury. Notify the dental depot or the manufacturer.

Potential explosion hazard

Do **not** use this product in areas subject to explosion hazards.

Repair

Do **not** repair the ultrasonic handpieces and tips on your own.

The BA instrument is the result of a precise development and must **not** be disassembled by third parties. In the case of service or repair work, please send the instrument in an assembled state to B.A. International Ltd. or to a service center approved by B.A. International.

Ultrasonic tips, spare parts, and accessories

Use only original B.A. International parts or parts/tips approved by B.A. International. Safe operation is **not** guaranteed for parts/tips that have not been approved by B.A. International and the handpiece may be damaged.

Endontic work

Cardiac pacemakers and implanted defibrillators

Combination with treatment centers

Use a cofferdam in patients for endodontic work.

B.A. International recommends that ultrasonic handpieces are not to be used if the patient or user has a cardiac pacemaker and/or implanted defibrillator.

Use the B.A. International ultrasonic handpieces only in combination with a treatment center approved by B.A. International. There is a risk of injury when used in combination with treatment centers from other manufacturers.

If you have any questions, please contact your dental supplies distributor or the manufacturer.

3 Technical description

3.1 Task

The ultrasonic handpiece converts electrical energy into a translational movement and transfers it to the ultrasonic tip, thus oscillating the distal end of the ultrasonic tip.

3.2 Ultrasonic handpiece design

Structure BACS125 / BACS126

A	Tip mount
B	Light aperture (only for BACS126)
C	Control ring
D	Hose nut

3.3 Technical Data

Ultrasonic handpiece

	BACS125	BACS126
Oscillation frequency in kHz	25 - 32	25 - 32
Spray function	x	x
Spray water regulation at the coupling	x	x
Light function	-	x
Amplitude of the instrument tip in μm ($100\mu\text{m} = 0.10\text{mm}$)	~ 180	~ 180

English

Dental treatment center

		BACS125	BACS126 ¹
Lamp voltage in V, at the lamp		-	3.6 ±0.1
Lamp current in A			
	Halogen lamp	-	-
	LED	-	0.03
Spray water flow in ml/min (with fully-opened water regulator)		60	60
Spray water pressure in bar (without tip)		2 ±0.2	2 ±0.2
Recommended water content in spray in ml/min		> 30 ml/min	> 30 ml/min

¹ depending on the treatment center and its accessories

Operating conditions

Ambient temperature	10 °C - 40 °C (50 °F - 104 °F)
Relative humidity	30 % - 95 %
Air pressure	700 hPa - 1060 hPa

This product bears the CE mark in accordance with the provisions of the Council Directive 93/42/EEC of June 14, 1993 concerning medical devices (MDD).

English

4 Pre-reprocessing

4.1 Initial start-up and longer breaks in use

- An adequate normalization phase must be maintained following transport before initial operation.
- Sterilize the ultrasonic handpiece and accessories prior to startup [→ 37].
- Clean and maintain the ultrasonic handpiece after longer breaks in use.

4.2 Prior to starting the work day

- Flush the water paths for 2 minutes.

4.3 Before each patient

1. Flush the water paths for 30 seconds.
2. Attach the ultrasonic handpiece [→ 21].

3. Insert the ultrasonic tip [→ 23].
 4. Check the ultrasonic tip for usability.
 5. Set a sufficient amount of cooling water (> 30 ml/min) [→ 25].
 6. Use filtered water only (< 50 µm).
 7. Check the nozzles for blockages and lime deposits, for example, and clean the nozzles if necessary.
 8. Make sure the ultrasonic tip fits tightly.
 9. Perform an axial load test of the ultrasonic tip.
- Tip:** Use a cofferdam in patients for endodontic work.

5 Operation

NOTE! Tips wear out and can bend resulting in a performance loss and risk of breakage. Use only ultrasonic tips that are **not** worn or bent. Do **not** reshape bent tips.

⚠ CAUTION! A loose ultrasonic tip can come loose from the thread or break off. This may cause injury! Use the ultrasonic handpiece only if the ultrasonic tip is firmly clamped in place.

⚠ CAUTION! Insufficient cooling leads to overheating of the preparation site, the tip, and the handpiece. Ensure that the water flow rate is > 30 ml/min.


⚠ CAUTION! You can puncture your hand with an ultrasonic handpiece when the tip is inserted. This may cause injury! When the ultrasonic handpiece is not in use, remove the tip or attach the torque wrench for protection.

⚠ CAUTION! Potentially hazardous optical radiation may cause harm to the eyes. Do **not** stare at the LED for longer periods of time while in operation.

5.1 Connecting the ultrasonic handpiece to the supply hose

Attaching ultrasonic handpiece

1. Compare the hole pattern in the handpiece with the small tubes in the hose coupling.
2. Hold the hose nut (A) firmly and insert the ultrasonic handpiece until it clicks audibly into place.


5.2 Detaching the ultrasonic handpiece from the supply hose

⚠ CAUTION! Do **not** detach the ultrasonic handpiece from the supply hose while it is in operation. This may cause injury!

Detaching ultrasonic handpiece

- ✓ The ultrasonic handpiece is **not** in operation.
- 1. Hold the supply hose by the control ring (B) and **not** by the hose nut (A)!
- 2. Detach the ultrasonic handpiece from the supply hose. Do **not** pull on or hold the supply hose.


5.3 Inserting and removing the ultrasonic tip

Inserting the ultrasonic tip

- ✓ The ultrasonic handpiece is **not** in operation.
- 1. Insert the ultrasonic tip into the torque wrench (A).
- 2. Slowly screw the ultrasonic tip into the torque wrench counterclockwise. Turn a quarter revolution past the resistance.
- 3. Detach the torque wrench.
- 4. Check the ultrasonic tip.

⚠ CAUTION! A loose ultrasonic tip can come loose from the thread or break off. This may cause injury! Use the ultrasonic handpiece only if the ultrasonic tip is firmly clamped in place.

Removing the ultrasonic tip


- ✓ The ultrasonic handpiece is **not** in operation.
- 1. Attach the torque wrench to the ultrasonic tip.
- 2. Unscrew the ultrasonic tip by turning it clockwise.

5.4 Adjusting the cooling spray

- Adjust the flow rate of the cooling water using the control ring (A) (> 30 ml/min).

Use the control ring (A) to regulate the water quantity in the BACS125 / BACS126. Turn the ring toward the + sign until perceptible resistance indicates maximum water flow is set. Turning to the - sign reduces the water flow until it is completely shut off.


Tip: You can measure the amount of cooling water with a measuring cup and watch.

5.5 Treatment

⚠ CAUTION! When the ultrasonic handpiece is not in use, detach the ultrasonic tip and attach the torque wrench to protect the tip. Make sure to keep the instrument away from soft tissues such as cheeks, lips, tongue, etc. during treatment. Contact causes friction that can lead to burns!

Select the most suitable instrument tip for each treatment.


IMPORTANT: Change tips only when the device is switched off.
[→ 23]

The illustration to the left shows how the handpiece should be held during treatment.

Always keep the instrument in motion, either up and down the tooth or across the proximal surface in lingual or buccal direction away from the tooth.

⚠ CAUTION! Always hold the instrument so that the motion of the tip is parallel to the tooth surface. Only slight pressure should be applied.

English

NOTE! During treatment, do not allow the device to come into contact with metal/ceramic restoration. Prostheses may otherwise become loose due to the high-frequency ultrasonic vibrations.

Since the instrument causes only slight pain, it is also ideal for the treatment of acute, painful periodontitis.

Use effective suction to ensure good visibility of the treatment area.

6 Post-reprocessing

6.1 After each treatment session

NOTE! Condition immediately, or at the latest, one hour after treatment.

⚠ CAUTION! You can puncture your hand with an ultrasonic handpiece when the tip is inserted. This may cause injury! When the ultrasonic handpiece is not in use, remove the tip or attach the torque wrench for protection.

- ✓ The ultrasonic handpiece is not in operation.
- ✓ Wear appropriate protective clothing.
- 1. Purge the water paths directly at the treatment unit for 30 seconds.
- 2. Predisinfect directly at the treatment center [→ 31].
- 3. Unscrew the ultrasonic tip using the torque wrench.

4. Detach the ultrasonic handpiece.
5. Transport the ultrasonic handpiece, ultrasonic tip and torque wrench to the hygiene room in a suitable transport container.
6. Conduct automatic conditioning. Manual conditioning [→ 34] is possible in exceptional cases if the national/local regulations are followed.
7. Sterilize the ultrasonic handpiece and accessories [→ 37].


7 Reprocessing

7.1 Conduct pre-disinfection

NOTE! Do not use any strong acidic or alkaline solutions ($5 < \text{pH} < 9$) or solutions containing chlorine.

- ✓ Wear appropriate protective clothing.
- ✓ All disinfectants must be approved in your country and have proven bactericidal, fungicidal and virucidal properties. Use only disinfectants with **no** protein-fixing properties.

English

1. Spray the surface with disinfectant.
 2. Wipe the disinfectant off with a cloth.
- ⚠ For further reprocessing, the instrument should be dry and free of residue.

Please observe the manufacturer's instructions for using instrument disinfectants.

7.2 Automated cleaning and disinfection

NOTE! Do not use any strong acidic or alkaline solutions ($5 < \text{pH} < 9$) or solutions containing chlorine.

IMPORTANT: Unscrew instrument tip and blow the remaining water out of the handpiece and tip.

The BACS125 / BACS126, torque wrench and ultrasonic tips can also be cleaned and disinfected in suitable cleaning and disinfection equipment.

NOTE! The ultrasonic tips can be cleaned and disinfected on the outside **only after removing them** from the torque wrench.

The cleaning and disinfection equipment used must be approved by its manufacturer for the cleaning and disinfection of dental instruments and comply with EN ISO 15883-1 (e.g. 95°C (203°F) and 10 min. holding time).

For further details, refer to the operating instructions supplied with the unit.

- ✓ The instrument is conditioned with a cleaning and disinfection device.
- 1. Check whether the instrument is clean after reprocessing under good lighting (min. 500 lux) and color rendering index (min. 80 Ra).

2. If they are still dirty, repeat the process.
 - ↳ For further reprocessing, the instrument should be dry and free of residue.
3. Blow the instrument out with max. 3 bar.
4. Pack the instrument in packaging material suitable for sterilization and storage. e.g. paper/plastic composite packaging.
5. Conduct sterilization [→ 37].

7.3 Manual cleaning and disinfection

IMPORTANT: Manual conditioning is possible in exceptional cases if the national/regional regulations are followed. The national/regional regulations are to be checked before.

NOTE! Condition immediately, or at the latest, one hour after treatment.

NOTE! Never clean in an ultrasound bath!


NOTE! Never immerse in disinfectant solution!

NOTE! Do not use any strong acidic or alkaline solutions ($5 < \text{pH} < 9$) or solutions containing chlorine.


IMPORTANT: Unscrew instrument tip and blow the remaining water out of the handpiece and tip.

IMPORTANT: Use a soft, clean, and disinfected brush for cleaning.

- ✓ Wear appropriate protective clothing.
- ✓ All disinfectants must be approved in your country and have proven bactericidal, fungicidal and virucidal properties. Use only disinfectants with **no** protein-fixing properties.

1. Brush the ultrasonic handpiece, detached tip and torque wrench under running water ($< 38\text{ }^{\circ}\text{C}$, $< 100\text{ }^{\circ}\text{F}$, at least drinking water quality) and good lighting (min. 500 lux) and color rendering index (min. 80 Ra) until no more dirt can be seen, for at least 10 seconds.
2. Blow spray channels with 2.5 - 3 bar until no more moisture exits, but at least 10 seconds.
3. Conduct thermal disinfection or unwrapped steam sterilization.
4. Pack the instruments in packaging material suitable for sterilization and storage. e.g. paper/plastic composite packaging.
5. Conduct sterilization [→ 37].

7.4 Sterilizing

- ✓ The ultrasonic handpiece is cleaned and disinfected.
- ✓ If necessary, the ultrasonic handpiece and accessories can be sterilized in packaging suitable for sterilization and storage, e.g. paper/plastic composite packaging or container.
- Sterilize the ultrasonic handpiece and accessories in the steam sterilizer with saturated water vapor.

English

Overpressure:	2.04 bar (29.59 psi)
Temperature:	134 °C (274 °F)
Holding time:	3 min.

Steam sterilizers that comply with the requirements of either EN 13060, class B or EN 13060, class S and are also suitable for the sterilization of ultrasonic handpieces and ultrasonic tips are approved.

NOTE! Do not exceed 140 °C (284 °F), even during the drying phase.

After sterilizing

1. Remove the instruments from the steam sterilizer immediately.
⚠ CAUTION! The instruments are hot. Risk of burns!
NOTE! Do **not** attempt to accelerate the cooling process by immersing the instruments in cold water. This can damage your instruments!

- 2.** Store all instruments so that they are protected from contamination.
- 3.** Sterilize the instruments again once the storage period has elapsed.


8 Maintenance


NOTE! Never oil the ultrasonic handpiece!

8.1 Maintaining the water channel

⚠ CAUTION! Insufficient cooling leads to overheating of the preparation site, the tip, and the handpiece. Ensure that the water flow rate is > 30 ml/min.

1. Use the cleaning wire to carefully clear the water channel of the tip regularly.
2. Blow the water channel out with max. 3 bar.

8.2 Cleaning the light guide surfaces


1. In order to avoid scratching the surfaces, blow off any dirt particles with air from a dry syringe.
2. Wipe the surfaces with a Q-tip or a soft cloth and alcohol.

8.3 Replacing the spray water cartridge

Required accessories

- Spray water cartridge
- B.A. Ultimate Spray

If the spray water cartridge leaks or is clogged.

**... at BACS125 / BACS126**


1. Detach the ultrasonic handpiece from the supply hose.
2. Detach the control ring from the hose coupling.
3. Remove the spray water cartridge (A).
4. Spray oil lightly onto the new spray water cartridge.
5. Insert the new spray water cartridge up to the stop.
6. Attach the control ring with the smooth ring side first onto the supply hose.

8.4 Replacing the halogen lamp/LED

... at BACS126

⚠ CAUTION! The lamp can be hot. Risk of burns! Allow the lamp to cool down.

1. Detach the ultrasonic handpiece.
2. Detach the control ring (B).
3. Insert a probe into the hole behind the lamp socket and push the defective halogen/LED lamp (A) out of the socket.
4. Insert the new halogen lamp (dark brown base/LED). Ensure the proper position of the contact surfaces.
NOTE! Pressing on the lens can destroy the LED. Use the assembly tool (C) to insert the LED.
5. Wipe the bulb of the halogen lamp with a clean cloth.
6. Attach the control ring with the smooth ring side first onto the supply hose.

If the LED does not light up:

- Remove the LED and re-insert it rotated 180° around the same axis.

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