

 HENRY SCHEIN®

 ZAHN™ LABORATORY  
DIVISION

# GIBRALTAR-ARTICULATION

## Plaster for Articulation, Type I

Escayola de precisión para articuladores, Tipo I

Plâtre pour Mise en Articulation, Type I

HS-Artikulationsgips, Typ I

Gesso di precisione per articolatori, Classe I

- 900-2876 Natural/Natural/Naturel/Natur/Naturale
- 900-2877 Synthetic/Sintética/Synthetiques/Synthetisch/Sintetico

For Professional Use Only  
Unicamente Para Uso Profesional  
Pour usage professionnel uniquement  
Nur für den Dentalgebrauch  
Solo per uso professionale

Contains:  
Contiene:  
Contient:  
Inhalt :  
Contenuto:

# 20kg

### Technical data

|  | Natural    | Synthetic  |
|--|------------|------------|
| Setting time/Tiempo de fraguado/Temps de prise/Gesamtaushärtezeit/Tempo di presa   | 4-5min     | 3-5min     |
| Compressive strength 1hr/Fuerza compresiva 1h/Resistance de compression 1 hour/Druckfestigkeit nach 1h/Resistenza alla compressione dopo 1 ora | 20         | 20         |
| Expansion 2 hrs/Expansión 2 horas/Expansion 2 heures/Expansion nach 2h/Espansioni di presa dopo 2 ore  | 0.04%      | 0.04%      |
| Mix ratio/Proporción de mezcla/Proportion de mélange/Mischungsverhältnis/Miscelatura   | 40ml/ 100g | 30ml/ 100g |

Distributed by:  
Distribuí par:

Distribuido por:  
Vertrieb durch:  
Distribuito da:

MADE IN GERMANY  
Fabricado en Alemania/Fabriqu  en Allemagne  
In Deutschland hergestellt/Prodotto in Germania

**HENRY SCHEIN INC.**  
Melville, NY 11747 USA  
Henry Schein U.K. Holdings Ltd.  
Gillingham ME8 0SB U.K.